


DURAN® SCREW CAP GL 45

with

2 HOSE CONNECTORS FOR
LABORATORY GLASS BOTTLES




DURAN
WHEATON
KIMBLE

Excellence in your hands

DURAN® SCREW CAP GL 45 with 2 HOSE CONNECTORS

The cap connection system has been specially developed for the easy connection of flexible tubing to the DURAN® GL 45 bottles.

The two polypropylene through connectors have angled top and straight underside connectors. The cap is ideal for use with soft elastic tubing that has an inner diameter of 6–9 mm (Viton®, Tygon®, silicone). An optional venting connector is available with, or without a syringe filter.



2-port GL 45 cap (29 310 28 07)

① Female Luer slip

② Syringe Filter Connector (11 298 29)

③ 5.8 mm connector with o-ring seal

④ Connector suitable for tubing with 6–9 mm inner diameter

⑤ 5.8 mm inner diameter

⑥ Connectors suitable for tubing with 6–9 mm inner diameter

⑦ Syringe Filter with connector (11 298 25)

Applications

- Possible biotech applications include the transfer of sterile media from one container to another using a peristaltic pump.

Benefits

- Highly versatile as the screw cap is based on the standard GL 45 thread
- Grey polypropylene cap center rotates freely, allowing the bottle to be unscrewed without the need to disconnect the tubing
- The angled top connectors prevent kinks in the tubing
- Temperature resistant up to +140 °C
- Fully autoclavable / dishwasher-proof



Ordering information

Description	Pack Unit	Cat. No.
DURAN® Screw Cap GL 45 with 2 hose connections EDPM Gasket	2	29 310 28 07
Syringe Filter 0.2 µm with connector	2	11 298 25
Non-sterile 0.2 µm filter and female Luer Slip to 5.8 mm male connector with o-ring seal		
Syringe Filter Connector only (without syringe filter)	2	11 298 29
Female Luer Slip to 5.8 mm male connector with o-ring seal		

You can obtain these products from your specialist laboratory ware dealer.



DWK Life Sciences GmbH
Hattenbergstrasse 10
55122 Mainz
Germany

Phone: +49 6131 - 1445 4131
Fax: +49 6131 - 1445 4016
sales@DWK-LifeSciences.com
www.DWK-LifeSciences.com

